# [www.sluhcentr.ru](http://www.sluhcentr.ru)

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- 1 On-Off Switch (travel lock)
- 2 Left-Both-Right-Switch
- 3 Left (Blue) Indicator Light
- 4 Right (Red) Indicator Light
- 5 Comfort-Clarity Balance\*
- 6 Volume Control
- 7 Home Button
- 8 Telephone Button
- 9 Program Change Button
- 10 LearnNow™ button inactive with Next™.
- \* Available in Next 16 only

## Legend

Creating innovative hearing solutions that make everyday listening easier for people with hearing loss.

[www.unitronhearing.com](http://www.unitronhearing.com)

## Smart Control Quick Reference

**next**™

Congratulations on purchasing your new Next™ hearing system with Smart Control. Smart Control enables you to personalize your listening experience in every situation so it's just right for you. See your Smart Control User Guide for step-by-step usage instructions.

### Here are a few tips to help you personalize your Next hearing system.

Use the + and - buttons to adjust your volume.

Use the 🏠 button to reset adjustments you made and set the HI back to the start program.

Use the ↺ program toggle to switch listening programs manually.

Use the 📞 button to directly access the designated program.

## Personalize in Two Ways

What Do You Want to Do?	Change the Quality of Sounds in a Noisy Situation	Change the Overall Sound Level
Which Control to Adjust:	<b>Comfort-Clarity Balance*</b>	<b>Volume Control</b>
How to use it:	<p>↗ If you want to sharpen speech, turn up to improve clarity</p> <p>↘ If you want to make a noisy environment comfortable, turn down to improve comfort</p>	<p><b>+</b> If the overall sound level is too soft, press the + button to increase loudness</p> <p><b>-</b> If the overall sound level is too loud, press the - button to decrease loudness</p>

\* Available in Next 16 only